

Pan Med Jr

Product Information

DESCRIPTION

As the moisture either drips on or flows over the **Pan Med Jr** geltab, it forms a gel, locking the geltab into place. This locking action never allows the geltab to move or block the drain line. **Pan Med Jr** geltabs uniquely designed formula with AT-8, turns the geltab into a very slow releasing gel, giving it a positive, long lasting and aggressive cleaning action. The geltab will clean sludge and other accumulations in condensate drain pans and lines, eliminating plugged condensate drain lines and stale water odors.

FEATURES AND BENEFITS

- Locks into place – never blocking a drain.
- Slow and continuous release - months of cleaning action
- Works till completely used up – nothing to remove or throw away.
- Eliminates odors from A/C – Stops foul and musty odors.

Available Sizes:

5112012- 1.5 Ton **Pan Med Jr**.

5112034- 3 Ton **Pan Med Jr**.

5112056- 5 Ton **Pan Med Jr**.

APPLICATIONS

For use in A/C condensate pans

PHYSICAL CHARACTERISTICS

Appearance	Amorphous
Crystalline	Solid
Odor	Aromatic
pH	7.4 ± .3 @ 10% Aqueous Solution
Specific Gravity	0.981 g/ml
Flash Point (°F)	Above 255 °F
Boiling Point	Not Applicable
Bio-Degradable	No

MARKETS

Parking Garages	Body Shops
Fleet Maintenance	Shopping Malls
Office Parks	Apartments
Restaurants	Heavy Equipment
Hospitals	Casinos
Nursing Homes	Prisons
Universities/Schools	And More!

DIRECTIONS

Place one **Pan Med Jr** geltab into the condensate drain pan and position so moisture drips on or runs over it, locking the geltab into place and slowly releasing it's chemical formula for a positive prolonged cleaning action.

**HARPER
SUPPLY**



7924 Camp Bowie West Blvd. • Fort Worth, TX 76116
817.529.1091 Office • 817.529.1092 Fax
www.HarperSupply.com

FOR PROFESSIONAL USE ONLY