

# Total-Kleen

Product Information

1842112

## DESCRIPTION

**Total-Kleen** is a powerful aerosol that thoroughly eliminates many harmful viruses and common microorganisms from hard, non-porous surfaces and leaves them clean and disinfected.

**Total-Kleen** is great for many places not just hospitals, but it will give you that hospital **Total-Kleen** effect.

## FEATURES AND BENEFITS

- A broad-spectrum disinfectant eliminates bacteria, viruses, molds and mildew.
- Fast and easy to use.
- Quickly disinfects hard, non-porous surfaces.
- Effective against: *Aspergillus niger*, HIV-1, Polio Virus type 1 and 2, Tubercule bacilli, and much more!

## PHYSICAL CHARACTERISTICS

Appearance	Aerosol
Odor	N/D
pH	N/D
Freezing Point	N/D
Flash Point	-156 °F
Specific Gravity/Density	0.79 g/ml estimated
Vapor Pressure	6.1 – 6.8 atm @ 70°F

## MARKETS

Medical Facilities	Nursing Homes
Government	Buildings
Cruise Ships	Air Ports
Hotels	Resorts
Schools	Industrial Areas

## DIRECTIONS

**To Disinfect:** After surface is cleaned, hold can upright 6-8 inches from surface to be treated. Spray 3-4 seconds until surface is wet. Allow to dry without wiping or allow spray to contact treated surface for a minimum of ten minutes prior to wiping.

**To Deodorize:** Hold can upright and point opening away from face. Spray into the air for 3-5 seconds directing spray toward center of room.

**To Control Surface Mildew:** After surface is cleaned, cover surface with mist. Repeat if needed.

**HARPER  
SUPPLY**



7924 Camp Bowie West Blvd. • Fort Worth, TX 76116  
817.529.1091 Office • 817.529.1092 Fax  
[www.HarperSupply.com](http://www.HarperSupply.com)

*FOR PROFESSIONAL USE ONLY*